

NEWS



FOR IMMEDIATE RELEASE

From: Jim Robbins
jim.robbins@brown-machine.com

For: Brown Machine LLC
330 N. Ross Street
Beaverton, MI 48612

877.702.4142 Toll Free
989.435.2190 Fax

Contact: Jim Robbins
989.435.2190 Ext.1216
989.859.7426 Cell

BROWN LR-2000 AND LR-2000S LIP ROLLERS: WHERE PERFORMANCE AND TECHNOLOGY CONVERGE

Beaverton, MI -- The new LR-2000S Lip Roller is designed to process containers from such advanced materials as polypropylene. Based on the highly successful 4-screw LR-2000 Lip Roller, this unit has successfully rolled 7-ounce polypropylene drink cups at over 2,500 cups per minute.

Brown's Model LR-2000 and LR-2000S Lip Rollers are designed to roll or curl lips on round, thermoplastic containers ranging in size from 1.5" (38 mm) to 5.5" (140 mm). The product is fed onto a conveyor belt which delivers the product to the crowd rollers. The crowd rollers spin the stack of parts and feeds them through the oven section to the curling screws. The curling screws capture the parts, roll the lip, and feed them down out of the machine.

Pre-formed containers can be manually or automatically fed to the lip roller by Brown's customized transfer conveyor. This system can process several thermoplastic containers that share a common stack height and top lip configuration, even though they may have different container diameters, and be lip-rolled with a common set of curling screws.

NEWS



FOR IMMEDIATE RELEASE

From: Jim Robbins
jim.robbsins@brown-machine.com

For: Brown Machine LLC
330 N. Ross Street
Beaverton, MI 48612

877.702.4142 Toll Free
989.435.2190 Fax

Contact: Jim Robbins
989.435.2190 Ext.1216
989.859.7426 Cell

Brown lip rollers are available as stand-alone units, or incorporated into a fully integrated thermoforming / lip-rolling system.

With maximum speed of 3,500 containers lip-rolled per minute, and product changeover time of two minutes or less, the Brown LR-2000 Series Lip Roller sets the industry standard for productivity and quality.

###